

Package ‘LipidMS’

January 19, 2021

Type Package

Title Lipid Annotation for LC-MS/MS DDA or DIA Data

Version 2.0.0

Author M Isabel Alcoriza-Balaguer, J Carlos Garcia-Canaveras

Maintainer M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Description Lipid annotation in untargeted liquid chromatography mass spectrometry lipidomics based on fragmentation rules. Alcoriza-Balaguer MI, Garcia-Canaveras JC, Lopez A, Conde I, Juan O, Carretero J, Lahoza A (2019) <doi:10.1021/acs.analchem.8b03409>.

Depends R (>= 4.0), shiny, utils

License GPL (>= 2)

LazyData TRUE

RoxygenNote 7.1.1

Imports readMzXmlData, CHNOSZ, scales, shinythemes, stats, graphics, grDevices

Suggests knitr, rmarkdown

VignetteBuilder knitr

Encoding UTF-8

NeedsCompilation yes

Repository CRAN

Date/Publication 2021-01-19 15:30:11 UTC

R topics documented:

adductsTable	3
assignDB	4
baconjdb	4
badb	5
carnitinesdb	5
CEdb	6

cerdb	6
chainFrgs	7
checkClass	9
checkIntensityRules	10
cldb	12
coelutingFrgs	13
coelutionScore	14
combineChains	14
confLevels	16
createLipidDB	16
crossTables	17
dataProcessing	18
ddaFrgs	20
dgdb	21
fadb	21
fahfadb	22
findCandidates	22
getInclusionList	24
hfadb	25
idBAneg	25
idCarpos	27
idCEpos	29
idCerneg	31
idCerpos	34
idCLneg	36
idDGpos	38
idFAHFAneg	41
idFAneg	43
idLPCneg	45
idLPCpos	47
idLPEneg	49
idLPEpos	51
idLPGneg	53
idLPIneg	55
idLPSneg	57
idMGpos	59
idNEG	61
idPCneg	62
idPCpos	65
idPENeg	67
idPEpos	70
idPGneg	72
idPGpos	74
idPIneg	77
idPOS	79
idPSneg	81
idSMpos	83
idSphneg	85

idSphPneg	87
idSphpos	89
idSphPpos	91
idTGpos	93
LipidMSapp	96
lysopadb	97
lysopcdb	97
lysopedb	98
lysopgdb	98
lysopidb	99
lysopsdb	99
mgdb	100
nlsphdb	100
organizeResults	101
padb	102
pcdb	103
pedb	103
pgdb	104
pidb	104
plotLipids	105
psdb	106
searchIsotopes	106
smdb	107
sphdb	107
sphPdb	108
tgdb	108
Index	109

adductsTable	<i>Adducts table</i>
--------------	----------------------

Description

Table of possible adducts to be employed by LipidMS and related information.

Usage

```
data("adductsTable")
```

Format

Data frame with 18 observations and the following 4 variables.

adduct character vector with the adducts names.

mdiff numeric vector indicating the mass differences.

charge numeric vector indicating the charge.

n numeric vector. It indicates if the ion is a monomer (1), a dimer (2), etc.

assignDB	<i>Load LipidMS default data bases</i>
----------	--

Description

load all LipidMS default data bases required to run identification functions.

Usage

```
assignDB()
```

Value

list of data frames

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:  
dbs <- assignDB()  
  
## End(Not run)
```

baconjdb	<i>Bile acids conjugates database</i>
----------	---------------------------------------

Description

Common bile acids conjugates. It can be modified to look for other BA species.

Usage

```
data("baconjdb")
```

Format

Data frame with 2 observations and the following 2 variables.

total character vector indicating the names of the conjugates.

Mass numeric vector with the neutral masses of the conjugates fragments.

badb*Bile acids database*

Description

In silico generated database for common bile acids.

Usage

```
data("badb")
```

Format

Data frame with 9 observations and the following 5 variables.

formula character vector with the molecular formulas.

total character vector containing the names of the BAs (i.e. CA, TDCA, GLCA...).

Mass numeric vector with the neutral masses.

conjugate character vector containing the conjugate of each BA.

base character vector containing the core of each BA.

carnitinesdb*Carnitines database*

Description

In silico generated database for common carnitines.

Usage

```
data("carnitinesdb")
```

Format

Data frame with 30 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

CEdb

CEs database

Description

In silico generated database for common CEs.

Usage

```
data("CEdb")
```

Format

Data frame with 30 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

cerdb

Ceramides database

Description

In silico generated database for common ceramides.

Usage

```
data("cerdb")
```

Format

Data frame with 52 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

chainFragments	<i>Search of chain specific fragments</i>
----------------	---

Description

Search of specific fragments that inform about the chains structure.

Usage

```
chainFragments(coelfragments, chainfragments, ppm = 10, candidates, f = NULL, dbs)
```

Arguments

coelfragments	coeluting fragments for each candidate. Output of coelutingFragments .
chainfragments	character vector containing the fragmentation rules for the chain fragments. If it is an empty vector, chains will be calculated based on the difference between the precursor and the other chain. See details.
ppm	m/z tolerance in ppm.
candidates	candidates data frame. If any chain needs to be calculated based on the difference between the precursor and the other chain, this argument will be required. Output of chainFragments .
f	known chains. If any chain needs to be calculated based on the difference between the precursor and the other chain, this argument will be required. Output of chainFragments .
dbs	list of data bases required for the annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be changed. If data bases have been customized using createLipidDB , they also have to be modified here.

Details

The chainfragments argument must contain the fragmentation rules which inform about the chains structure. For example, in the case of PG subclass, the chain in sn1 position is identified by the lysoPG as M-H resulting from the loss of the FA chain of sn2; and the chain in sn2 position is identified as the free FA chain as M-H. These two fragments need to be searched in two different steps: in the first step we will look for lysoPGs coeluting with the precursor using chainfragments = c("lysopg_M-H"); then, we will look for FA chains using chainfragments = c("fa_M-H"). This information can be combined later using [combineChains](#) function.

To indicate the fragments to be searched, the class of lipid is written using the same names as the LipidMS databases without the "db" at the end (i.e. pa, dg, lysopa, mg, CE, etc.), and the adduct has to be indicated as it appears in the adductsTable, both parts separated by "_". In case some chain needs to be searched based on a neutral loss, this can be defined using "NL-" prefix, followed by the database and adduct. If this neutral loss is employed to find the remaining chain, "cbdiff-" prefix allows to calculate the difference in carbons and doubles bounds between the precursor and the building block found. For example, "cbdiff-dg_M+H-H2O" will look for DG as M+H-H2O and

then, it will return the difference between their number of carbons and double bounds and the ones from the precursor. Otherwise, "NL-mg_M+H-H2O" will look for fragments coming from the loss of MGs.

In case these fragments identified as losses from the precursors are going to be employed for the intensity rules, this same prefix has to be added.

If a chain is calculated based on the difference of total number of carbons and double bounds between the precursor and a previously searched chain, chainfrags argument must be a character vector c("") and candidates data frame and chain fragments list must be provided.

Value

List of data frames with the chain fragments found.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
dbs <- assignDB()

MS1 <- LipidMSdata2::msobjectDIAneg$peaklist$MS1
MS1 <- MS1[MS1$isotope %in% c("[M+0]"), !colnames(MS1) %in% c("isotope", "group")]
MS2 <- LipidMSdata2::msobjectDIAneg$peaklist$MS2[,c("m.z", "RT", "int", "peakID")]
rawData <- rbind(LipidMSdata2::msobjectDIAneg$MS1, LipidMSdata2::msobjectDIAneg$MS2)

candidates <- findCandidates(MS1 = MS1, db = dbs$pgdb, ppm = 10,
rt = c(0, 2000), adducts = c("M-H"), rttol = 10, dbs = dbs,
rawData = rawData, coelCutoff = 0.8)
coelfrags <- coelutingFrgs(candidates, MSMS, rttol = 10, rawData = rawData,
coelCutoff = 0.8)
classConf <- checkClass(candidates, coelfrags,
clfrags = c(227.0326, 209.022, 74.0359),
clrequisites = c(FALSE, FALSE, FALSE, FALSE),
ftype = c("F", "F", "NL"), ppm = 10, dbs = dbs)

sn1 <- chainFrgs(coelfrags, chainfrags = c("lysopg_M-H"),
candidates = candidates, ppm = 10, dbs = dbs)
sn2 <- chainFrgs(coelfrags, chainfrags = c("fa_M-H"), ppm = 10,
candidates = candidates, dbs = dbs)

## End(Not run)
```

checkClass	<i>Search of class fragments to confirm the lipid class.</i>
------------	--

Description

Search of characteristic fragments that confirm a given lipid class.

Usage

```
checkClass(candidates, coelfrags, clfrags, ftype, clrequisites, ppm = 10, dbs)
```

Arguments

candidates	output of findCandidates function.
coelfrags	list of peaks coeluting with each candidate. Output of coelutingFrgs .
clfrags	vector containing the expected fragments for a given lipid class. See details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See details.
clrequisites	logical vector indicating if each class fragment is required or not. If none of the fragment is required, at least one of them must be present within the coeluting fragments. If the presence of any fragment excludes the class, it can be specified by using "excluding".
ppm	m/z tolerance in ppm.
dbs	list of data bases required for the annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be changed. If data bases have been customized using createLipidDB , they also have to be modified here. It is employed when some fragment belongs to "BB" ftype.

Details

clfrags, ftype and clrequisites will indicate the rules to confirm a lipid class. All three arguments must have the same length.

This function allows three different types of fragments: fragments with a specific m/z as for example 227.0326 for PG in negative mode, which needs to be defined as clfrags = c(227.0326) and ftype = c("F"); neutral losses such as the head group of some PL (i.e. NL of 74.0359 in PG in negative mode), which will be defined as clfrags = c(74.0359) and ftype = c("NL"); or building blocks resulting from the loss of some groups, as for example, PA as M-H resulting from the loss of the head group (glycerol) in PG in ESI-, which will be defined as clfrags = c("pa_M-H") and ftype = c("BB"). The last two options could define the same fragments. In this case just one of them would be necessary.

When using the third type of fragment ("BB"), the building block will be specified in lower case (i.e. pa, dg, lysopa, mg, etc.) and the adduct will be given as it appears in the adductsTable, both separated by "_". Names for the building blocks are the ones used for the LipidMS databases without the "db" at the end.

In case the presence of a fragment indicates that the candidate does not belong to the lipid class (i.e. loss of CH₃ in PE, which corresponds to a PC actually), this will be specified by using `clrequisites = c("excluding")`.

Value

List with 2 elements: a matrix with logical values (presence/absence) of each expected fragment (columns) for each candidate (rows), and a logical vector with the confirmation of the lipid class for each candidate.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
dbs <- assignDB()

MS1 <- LipidMSdata2::msobjectDIAneg$peaklist$MS1
MS1 <- MS1[MS1$isotope %in% c("[M+0]"), !colnames(MS1) %in% c("isotope", "group")]
MS2 <- LipidMSdata2::msobjectDIAneg$peaklist$MS2[,c("m.z", "RT", "int", "peakID")]
rawData <- rbind(LipidMSdata2::msobjectDIAneg$MS1, LipidMSdata2::msobjectDIAneg$MS2)

candidates <- findCandidates(MS1 = MS1, db = dbs$pgdb, ppm = 10,
rt = c(0, 2000), adducts = c("M-H"), rttol = 10, dbs = dbs,
rawData = rawData, coelCutoff = 0.8)
coelfrags <- coelutingFragments(candidates, MSMS, rttol = 10, rawData = rawData,
coelCutoff = 0.8)

classConf <- checkClass(candidates, coelfrags,
clfrags = c(227.0326, 209.022, 74.0359),
clrequisites = c(FALSE, FALSE, FALSE, FALSE),
ftype = c("F", "F", "NL"), ppm = 10, dbs = dbs)

## End(Not run)
```

checkIntensityRules	<i>Check intensity rules</i>
---------------------	------------------------------

Description

Check intensity rules to confirm chains position.

Usage

```
checkIntensityRules(inrules, rates, intrequired, nchains, combinations)
```

Arguments

<code>inrules</code>	character vector specifying the fragments to compare. See details.
<code>rates</code>	character vector with the expected rates between fragments given as a string (i.e. "3/1"). See details.
<code>intrequired</code>	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
<code>nchains</code>	number of chains of the targeted lipid class.
<code>combinations</code>	output of combineChains .

Details

This function will be employed when the targeted lipid class has more than one chain.

Taking PG subclass as an example, intensities of lysoPG fragments (informative for sn1) can be employed to confirm the chains structure (`inrules = c("lysopg_sn1/lysopg_sn1")`). In this case, the intensity of the lysoPG resulting from the loss of the FA chain in sn2 is at least 3 times greater (`rates = c("3/1")`) than the lysoPG resulting from the loss of the FA chain in sn1.

For the `inrules` argument, "/" will be use to separate the fragments related to each chain (sn1/sn2/etc), and "_" will be use to indicate the list in which they'll be searched. This will depend on the chain fragments rules defined previously. Following the example, as we use lysoPG to define the sn1 position, both fragments will be searched in this list (sn1).

For classes with more than one FA chain, if some intensity rule should be employed to identify their position but they are no defined yet, use "Unknown". If it is not necessary because the fragmentation rules are informative enough to define the position (i.e. sphingolipid species), just leave an empty vector.

Value

List of logical vectors with the confirmation for each combination.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
dbs <- assignDB()

MS1 <- LipidMSdata2::msobjectDIAneg$peaklist$MS1
MS1 <- MS1[MS1$isotope %in% c("[M+0]"), !colnames(MS1) %in% c("isotope", "group")]
MS2 <- LipidMSdata2::msobjectDIAneg$peaklist$MS2[,c("m.z", "RT", "int", "peakID")]
```

```

rawData <- rbind(LipidMSdata2::msobjectDIAneg$MS1, LipidMSdata2::msobjectDIAneg$MS2)

candidates <- findCandidates(MS1 = MS1, db = dbs$pgdb, ppm = 10,
rt = c(0, 2000), adducts = c("M-H"), rttol = 10, dbs = dbs,
rawData = rawData, coelCutoff = 0.8)
coelfrags <- coelutingFragments(candidates, MSMS, rttol = 10, rawData = rawData,
coelCutoff = 0.8)
classConf <- checkClass(candidates, coelfrags,
clfrags = c(227.0326, 209.022, 74.0359),
clrequisites = c(FALSE, FALSE, FALSE, FALSE),
ftype = c("F", "F", "NL"), ppm = 10, dbs = dbs)
sn1 <- chainFragments(coelfrags, chainfrags = c("lysopg_M-H"),
candidates = candidates, ppm = 10, dbs = dbs)
sn2 <- chainFragments(coelfrags, chainfrags = c("fa_M-H"), ppm = 10,
candidates = candidates, dbs = dbs)
chainsComb <- combineChains(candidates, nchains=2, sn1, sn2)

intConf <- checkIntensityRules(inrules = c("lysopg_sn1/lysopg_sn1"),
rates = c("2/1"), intrequired = c(TRUE), nchains=2, chainsComb)

## End(Not run)

```

cldb

Cardiolipins database

Description

In silico generated database for commo CLs.

Usage

```
data("cldb")
```

Format

Data frame with 714 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

coelutingFrgs	<i>Coeluting fragments extraction</i>
---------------	---------------------------------------

Description

Given a RT and a list of peaks, this function subsets all coeluting fragments within a rt windows. It is used by identification functions to extract coeluting fragments from high energy functions for candidate precursor ions.

Usage

```
coelutingFrgs(
  precursors,
  products,
  rttol,
  rawData = data.frame(),
  coelCutoff = 0
)
```

Arguments

precursors	candidates data frame. Output of findCandidates .
products	peaklist for MS2 function (MSMS).
rttol	rt window in seconds.
rawData	raw scans data. Output of dataProcessing function (MSMS\$rawData).
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied.

Value

List of data frames with the coeluting fragments for each candidate.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
dbs <- assignDB()

MS1 <- LipidMSdata2::msobjectDIAneg$peaklist$MS1
MS1 <- MS1[MS1$isotope %in% c("[M+0]"), !colnames(MS1) %in% c("isotope", "group")]
MS2 <- LipidMSdata2::msobjectDIAneg$peaklist$MS2[,c("m.z", "RT", "int", "peakID")]
```

```

rawData <- rbind(LipidMSdata2::msobjectDIAneg$MS1, LipidMSdata2::msobjectDIAneg$MS2)

candidates <- findCandidates(MS1 = MS1, db = dbs$pgdb, ppm = 10,
rt = c(0, 2000), adducts = c("M-H"), rttol = 10, dbs = dbs,
rawData = rawData, coelCutoff = 0.8)

coelfrags <- coelutingFragments(candidates, MSMS, rttol = 10, rawData = rawData,
coelCutoff = 0.8)

## End(Not run)

```

coelutionScore	<i>calculate coelution score between two peaks</i>
----------------	--

Description

Calculate coelution score between two peaks.

Usage

```
coelutionScore(peak1, peak2, rawData)
```

Arguments

peak1	character vector specifying the peakID of the first peak.
peak2	character vector specifying the peakID of the second peak.
rawData	data frame with raw data for each scan. it need to have at least 5 columns: m.z, RT, int, Scan (ordinal number for a given MS function) and peakID (peakID to which it has been assigned). #’ @keywords internal

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

combineChains	<i>Combine chain fragments that could belong to the same precursor.</i>
---------------	---

Description

It calculates combinations of chain fragments that sum up the same number of carbons and double bounds as the precursor.

Usage

```
combineChains(candidates, nchains, sn1, sn2, sn3, sn4)
```

Arguments

candidates	candidates data frame. Output of findCandidates .
nchains	number of chains of the targeted lipid class.
sn1	list of chain fragments identified for sn1 position. Output of chainFragments .
sn2	list of chain fragments identified for sn2 position. Output of chainFragments . If required.
sn3	list of chain fragments identified for sn3 position. Output of chainFragments . If required.
sn4	list of chain fragments identified for sn4 position. Output of chainFragments . If required.

Value

List of data frames with candidate chains structures.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
dbs <- assignDB()

MS1 <- LipidMSdata2::msobjectDIAneg$peaklist$MS1
MS1 <- MS1[MS1$isotope %in% c("[M+0]"), !colnames(MS1) %in% c("isotope", "group")]
MS2 <- LipidMSdata2::msobjectDIAneg$peaklist$MS2[,c("m.z", "RT", "int", "peakID")]
rawData <- rbind(LipidMSdata2::msobjectDIAneg$MS1, LipidMSdata2::msobjectDIAneg$MS2)

candidates <- findCandidates(MS1 = MS1, db = dbs$pgdb, ppm = 10,
rt = c(0, 2000), adducts = c("M-H"), rttol = 10, dbs = dbs,
rawData = rawData, coelCutoff = 0.8)
coelfrags <- coelutingFragments(candidates, MSMS, rttol = 10, rawData = rawData,
coelCutoff = 0.8)
classConf <- checkClass(candidates, coelfrags,
clfrags = c(227.0326, 209.022, 74.0359),
clrequisites = c(FALSE, FALSE, FALSE, FALSE),
ftype = c("F", "F", "NL"), ppm = 10, dbs = dbs)
sn1 <- chainFragments(coelfrags, chainfrags = c("lysopg_M-H"),
candidates = candidates, ppm = 10, dbs = dbs)
sn2 <- chainFragments(coelfrags, chainfrags = c("fa_M-H"), ppm = 10,
candidates = candidates, dbs = dbs)

chainsComb <- combineChains(candidates, nchains=2, sn1, sn2)

## End(Not run)
```

confLevels	<i>Confidence Annotation Levels</i>
------------	-------------------------------------

Description

Confidence annotation levels and their hierarchy.

Usage

```
data("confLevels")
```

Format

Data frame with 5 observations and 2 variables.

level character vector with the names of the annotation levels.

order numeric vector that indicates the hierarchical order.

createLipidDB	<i>Customizable lipid DBs creator</i>
---------------	---------------------------------------

Description

It allows to create easy-customizable lipid DBs for annotation with LipidMS package.

Usage

```
createLipidDB(lipid, chains, chains2)
```

Arguments

lipid	character value indicating the class of lipid. See Details.
chains	character vector indicating the FA chains to be employed
chains2	character vector containing the sphingoid bases to be employed if required.

Details

lipidClass argument needs to be one of the following character values: "Cer", "CerP", "GlcCer", "SM", "Carnitine", "CE", "FA", "HFA", "Sph" (sphingoid bases), "SphP", "MG", "LPA", "LPC", "LPE", "LPG", "LPI", "LPS", "FAHFA", "DG", "PC", "PE", "PG", "PI", "PS", "PA", "TG", "CL" or "all".

Value

List with the requested dbs (data frames)

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafes.es>

Examples

```
fas <- c("8:0", "10:0", "12:0", "14:0", "14:1", "15:0", "16:0", "16:1",
"17:0", "18:0", "18:1", "18:2", "18:3", "18:4", "20:0", "20:1", "20:2",
"20:3", "20:4", "20:5", "22:0", "22:1", "22:2", "22:3", "22:4", "22:5",
"22:6", "24:0", "24:1", "26:0")
sph <- c("16:0", "16:1", "18:0", "18:1")
newdb <- createLipidDB(lipid = "PC", chains = fas, chains2 = sph)
```

crossTables

Cross the original MS1 peaklist with the annotation results

Description

Cross the original MS1 peaklist with the annotation results.

Usage

```
crossTables(MS1, results, ppm = 5, rttol = 10, dbs)
```

Arguments

MS1	data frame containing all peaks from the full MS function. It must have three columns: m.z, RT (in seconds) and int (intensity).
results	data frame. Output of identification functions.
ppm	mass tolerance in ppm.
rttol	rt tolerance to match peaks in seconds.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Value

Data frame with 6 columns: m.z, RT, int, LipidMS_id, adduct and confidence level for the annotation. When multiple IDs are proposed for the same feature, they are sorted based on the annotation level.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafes.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPCneg(LipidMSdata2::msobjectDIAneg)
crossTables(msobject$peaklist$MS1, msobject$results,
ppm = 10, rttol = 10)
## End(Not run)
```

dataProcessing

Process mzXML files: peakpicking and deisotoping

Description

Process mzXML files: peak-picking using enviPick and deisotoping using an adaptation of the CAMERA algorithm.

Usage

```
dataProcessing(
  file,
  acquisitionmode,
  polarity,
  dmzagglom = 10,
  drtagglom = 500,
  drtclust = 100,
  minpeak = c(5, 3),
  drtgap = 10,
  drtminpeak = 20,
  drtmaxpeak = 100,
  recurs = 5,
  sb = c(3, 2),
  sn = 2,
  minint = c(1000, 100),
  weight = c(2, 3),
  dmzIso = 10,
  drtIso = 5
)
```

Arguments

file	path of the mzXML input file.
acquisitionmode	character value: DIA or DDA (acquisition mode).
polarity	character value: negative or positive.

dmzagglom	m/z tolerance (in ppm) used for partitioning and clustering. 10 by default.
drtagglom	rt window used for partitioning (in seconds). 500 by default.
drtclust	rt window used for clustering (in seconds). 100 by default.
minpeak	minimum number of measurements required for a peak. By default, 5 for MS1 and 4 for MS2.
drtgap	maximum RT gap length to be filled (in seconds). 5 by default.
drtminpeak	minimum RT width of a peak (in seconds). 20 by default. At least minpeak within the drtminpeak window are required to define a peak.
drtmaxpeak	maximum RT width of a single peak (in seconds). 100 by default.
recurs	maximum number of peaks within one EIC. 5 by default.
sb	signal-to-base ratio. By default, 3 for MS1 and 2 for MS2.
sn	signal-to-noise ratio. 2 by default.
minint	minimum intensity of a peak. By default, 1000 for MS1 and 100 for MS2.
weight	weight for assigning measurements to a peak. Optional. By default, 2 for MS1 and 3 for MS2.
dmzIso	numeric. Mass tolerance for isotope matching. 10 by default.
drtIso	numeric. Time windows for isotope matching. 5 by default.

Details

This function executes 2 steps: 1) peak-picking based on `enviPick` package and 2) isotope annotation.

Value

a `msoject` that contains metadata of the `mzXML` file, raw data and extracted peaks.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@iislafe.es>

References

Peak-picking algorithm has been imported from `enviPick` R-package: <https://cran.r-project.org/web/packages/enviPick/index>.

Examples

```
## Not run:
msoject <- dataProcessing("input_file.mzXML", acquisitionmode="DIA", polarity,
dmzagglom = 25, drtagglom = 500, drtclust = 60, minpeak = c(5, 3),
drtgap = 5, drtminpeak = 20, drtmaxpeak = 100, recurs = 5, sb = c(3, 2),
sn = 2, minint = c(1000, 100), weight = 2, dmzIso = 10, drtIso = 5)

## End(Not run)
```

ddaFragments

*MS/MS scan extraction of a precursor in DDA***Description**

This function searches for the closest precursor selected for MS2 in DDA that matches m.z tolerance and RT window of a list of candidates and extracts their fragments.

Usage

```
ddaFragments(candidates, precursors, rawData, ppm)
```

Arguments

candidates	candidates data frame. Output of findCandidates .
precursors	data frame with the whole list of precursors selected for MS2.
rawData	peaklist for MS2 function (MSMS).
ppm	m/z tolerance in ppm.

Details

MS2 scans for a given precursor are searched within a rt window from minrt-rttol/2 to maxrt+rttol/2. If the same precursor was selected several times along the peak, the closest scan to the rt at the peak maximum is selected for annotation.

Value

List of data frames with the fragments for each candidate.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
dbs <- assignDB()

MS1 <- LipidMSdata2::msobjectDDAneg$peaklist$MS1
MS1 <- MS1[MS1$isotope %in% c("[M+0]"), !colnames(MS1) %in% c("isotope", "group")]
MS2 <- LipidMSdata2::msobjectDDAneg$peaklist$MS2[,c("m.z", "RT", "int", "peakID")]
rawData <- rbind(LipidMSdata2::msobjectDDAneg$MS1, LipidMSdata2::msobjectDDAneg$MS2)
precursors <- LipidMSdata2::msobjectDDAneg$metaData$scansMetadata[
  LipidMSdata2::msobjectDDAneg$metaData$scansMetadata$collisionEnergy > 0 &
  msobjectDDAneg$metaData$scansMetadata$msLevel == 2, c("retentionTime", "precursor", "Scan")]
```

```
candidates <- findCandidates(MS1 = MS1, db = dbs$cerdb, ppm = 10,
rt = c(0, 2000), adducts = c("M-H"),
rttol = 10, dbs = dbs, rawData = rawData_neg$rawScans, coelCutoff = 0.8)

coelfrags <- ddaFragments(candidates, precursors, rawData, ppm = 10)

## End(Not run)
```

dgdb

DGs database

Description

In silico generated database for common DGs.

Usage

```
data("dgdb")
```

Format

Data frame with 147 observations and the following 3 variables.

`formula` character vector containing molecular formulas.

`total` character vector indicating the total number of carbons and double bonds of the chains.

`Mass` numeric vector with the neutral masses.

fadb

FAs database

Description

In silico generated database for common FAs.

Usage

```
data("fadb")
```

Format

Data frame with 30 observations and the following 3 variables.

`formula` character vector containing molecular formulas.

`total` character vector indicating the total number of carbons and double bonds of the chains.

`Mass` numeric vector with the neutral masses.

fahfadb

FAHFAs database

Description

In silico generated database for common FAHFAs.

Usage

```
data("fahfadb")
```

Format

Data frame with 147 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

findCandidates

Search of lipid candidates of a certain class

Description

Search of lipid candidates from a peaklist based on a set of expected adducts.

Usage

```
findCandidates(  
  MS1,  
  db,  
  ppm,  
  rt,  
  adducts,  
  rttol = 3,  
  dbs,  
  rawData = data.frame(),  
  coelCutoff = 0  
)
```

Arguments

MS1	peaklist of the MS function. Data frame with 3 columns: m.z, RT (in seconds) and int (intensity).
db	database (i.e. pcdB, dgdb, etc.). Data frame with at least 2 columns: Mass (exact mass) and total (total number of carbons and double bound of the FA chains, i.e. "34:1").
ppm	m/z tolerance in ppm.
rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	character vector containing the expected adducts to search for (i.e. "M+H", "M+Na", "M-H", etc.). See details.
rttol	rt tolerance in seconds to match adducts.
dbs	list of data bases required for the annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be changed. If data bases have been customized using createLipidDB , they also have to be modified here.
rawData	raw scans data. Output of dataProcessing function (MS1\$rawData).
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied.

Details

[findCandidates](#) looks for matches between the m/z of the MS1 peaklist and the expected m/z of the candidates in the database for each adduct. If several adducts are expected, results are combined.

Adducts allowed are contained in adductsTable data frame, which can be modified if required (see [adductsTable](#)).

Value

Data frame with the found candidates. It contains 6 columns: m.z, RT, int (from the peaklist data.frame), ppms, cb (total number of carbons and double bounds of the FA chains) and adducts.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
dbs <- assignDB()

MS1 <- LipidMSdata2::msobjectDIAneg$peaklist$MS1
MS1 <- MS1[MS1$isotope %in% c("[M+0]"), !colnames(MS1) %in% c("isotope", "group")]
MS2 <- LipidMSdata2::msobjectDIAneg$peaklist$MS2[,c("m.z", "RT", "int", "peakID")]
```

```

rawData <- rbind(LipidMSdata2::msobjectDIAneg$MS1, LipidMSdata2::msobjectDIAneg$MS2)

candidates <- findCandidates(MS1 = MS1, db = dbs$pgdb, ppm = 10,
rt = c(0, 2000), adducts = c("M-H"), rttol = 10, dbs = dbs,
rawData = rawData, coelCutoff = 0.8)

# If any adduct is not in the adductsTable, it can be added:

adductsTable2 <- rbind(dbs$adductsTable,
data.frame(adduct = "M+X", mdiff = 44.9982, n = 1, charge = -1,
stringsAsFactors = FALSE))
dbs <- assignDB()
dbs$adductsTable <- adductsTable2

candidates <- findCandidates(MS1 = MS1, db = dbs$pgdb, ppm = 10,
rt = c(0, 2000), adducts = c("M-H", "M+X"), rttol = 10, dbs = dbs,
rawData = rawData, coelCutoff = 0.8)

## End(Not run)

```

getInclusionList

Obtain an inclusion list from the annotation results

Description

Obtain an inclusion list from the annotation results.

Usage

```
getInclusionList(results, adductsTable = LipidMS::adductsTable)
```

Arguments

results	data frame. Output of identification functions.
adductsTable	data frame with the adducts allowed and their mass difference.

Value

Data frame with 6 columns: formula, RT, neutral mass, m/z, adduct and the compound name.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPOS(LipidMSdata2::msobjectDIApos)
getInclusionList(msobject)

## End(Not run)
```

hfadb	<i>HFAs database</i>
-------	----------------------

Description

In silico generated database for common HFAs.

Usage

```
data("hfadb")
```

Format

Data frame with 30 observations and the following 3 variables.

`formula` character vector containing molecular formulas.

`total` character vector indicating the total number of carbons and double bounds of the chains.

`Mass` numeric vector with the neutral masses.

idBANeg	<i>Bile Acids (BA) annotation for ESI-</i>
---------	--

Description

BA identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idBAneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H"),
  conjfrag = c("baconj_M-H"),
  bafrag = c("ba_M-H-H2O", "ba_M-H-2H2O"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

msobject	an msobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursor and product ions. By default, 3 seconds.
rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for BA in ESI-. Adducts allowed can be modified in the adductsTable (dbs argument).
conjfrag	character vector containing the fragmentation rules for the BA-conjugates. By default just taurine and glycine are considered, but baconjdb can be modified to add more possible conjugates. See chainFrag s for details. It can also be an empty vector.
bafrag	character vector containing the fragmentation rules for other BA fragments. See chainFrag s for details. It can be an empty vector.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idBAneg function involves 3 steps. 1) FullMS-based identification of candidate BA as M-H. 2) Search of BA-conjugate fragments if required. 3) Search of fragments coming from the loss of H₂O.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (MS-only if no rules are defined, or Subclass level if they are supported by fragments) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idBAneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idCarpos

Acylcarnitine annotation for ESI+

Description

Acylcarnitines identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idCarpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H", "M+Na"),
  clfrags = c(60.0807, 85.0295, "fa_M+H-H2O"),
```

```

    clrequired = c(F, F, F),
    ftype = c("F", "F", "BB"),
    chainfrags_sn1 = c("fa_M+H-H2O"),
    coelCutoff = 0.8,
    dbs
  )

```

Arguments

<code>msoject</code>	an <code>msoject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for Carnitines in ESI+. Adducts allowed can be modified in <code>adductsTable</code> (<code>dbs</code> argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments. See chainFrgs for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>dbs</code>	list of data bases required for annotation. By default, <code>dbs</code> contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>dbs</code> may need to be supplied. See createLipidDB and assignDB .

Details

`idCarpos` function involves 3 steps. 1) FullMS-based identification of candidate carnitines as M+H and M+Na. 2) Search of carnitine class fragments: 60.0807 and 85.0295 or its loss (FA as M+H-H2O) coeluting with the precursor ion. 3) Search of specific fragments coming from the FA chain (FA as M+H-H2O).

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m/z error), confidenceLevel (in this case, as Carnitines only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idCarpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idCEpos

Cholesteryl Esters (CE) annotation for ESI+

Description

CE identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idCEpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("2M+NH4", "2M+Na", "M+NH4", "M+Na"),
  clfrags = c(369.3516, "fa_M+H-H2O"),
```

```

    clrequired = c(F, F),
    ftype = c("F", "BB"),
    chainfrags_sn1 = c("fa_M+H-H2O"),
    coelCutoff = 0.8,
    dbs
  )

```

Arguments

<code>mobject</code>	an <code>mobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for CE in ESI+. Adducts allowed can be modified in <code>adductsTable</code> (<code>dbs</code> argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments. See chainFrgs for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>dbs</code>	list of data bases required for annotation. By default, <code>dbs</code> contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>dbs</code> may need to be supplied. See createLipidDB and assignDB .

Details

`idCEpos` function involves 3 steps. 1) FullMS-based identification of candidate CE as $2M+NH_4$, $2M+Na$, $M+NH_4$ and $M+Na$. 2) Search of CE class fragments: 369.3516 or its loss (FA as $M+H-H_2O$) coeluting with the precursor ion. 3) Search of specific fragments that confirm chain composition (FA as $M+H-H_2O$).

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *mz*, RT (in seconds), *I* (intensity, which comes directly from de input), Adducts, ppm (*m.z* error), confidenceLevel (in this case, as CE only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idCEpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idCerneg

Ceramides (Cer) annotation for ESI-

Description

Cer identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idCerneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H", "M+CH3COO"),
  clfrags = c(),
```

```

    clrequired = c(),
    ftype = c(),
    chainfrags_sn1 = c("NL-nlsph_M-H", "sph_M-H-2H2O", "sph_M-H-H2O"),
    chainfrags_sn2 = c("fa_Mn-1.9918"),
    intrules = c(),
    rates = c(),
    intrequired = c(),
    coelCutoff = 0.8,
    dbs
  )

```

Arguments

<code>msobject</code>	an <code>msobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for Cer in ESI-. Adducts allowed can be modified in <code>adductsTable</code> (<code>dbs</code> argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFragments for details.
<code>chainfrags_sn2</code>	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFragments for details. If empty, it will be estimated based on the difference between precursors and sn1 chains.
<code>intrules</code>	character vector specifying the fragments to compare. See checkIntensityRules .
<code>rates</code>	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
<code>intrequired</code>	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>dbs</code>	list of data bases required for annotation. By default, <code>dbs</code> contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>dbs</code> may need to be supplied. See createLipidDB and assignDB .

Details

idCerneg function involves 5 steps. 1) FullMS-based identification of candidate Cer as M-H and M+CH₃COO. 2) Search of Cer class fragments: there are no class fragment by default. 3) Search of specific fragments that inform about the sphingoid base (Sph as M-H-2H₂O resulting from the loss of the FA chain or loss of part of the sphingoid base) and the FA chain (FA as M-H but with a N instead of an O, what means a mass difference of 1.9918 from the exact mass of the FA). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In this case, there are no intensity rules by default.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idCerneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idCerpos

*Ceramides (Cer) annotation for ESI+***Description**

Ceramides identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idCerpos(
  msubject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H-H2O", "M+Na", "M+H"),
  clfrags = c(),
  clrequired = c(),
  ftype = c(),
  chainfrags_sn1 = c("sph_M+H-2H2O"),
  chainfrags_sn2 = c(""),
  intrules = c(),
  rates = c(),
  intrequired = c(),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>msubject</code>	an <code>msobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for Cer in ESI+. Adducts allowed can be modified in <code>adductsTable</code> (<code>dbs</code> argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.

<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments in <code>sn1</code> position. See chainFrgs for details.
<code>chainfrags_sn2</code>	character vector containing the fragmentation rules for the chain fragments in <code>sn2</code> position. See chainFrgs for details. If empty, it will be estimated based on the difference between precursors and <code>sn1</code> chains.
<code>intrules</code>	character vector specifying the fragments to compare. See checkIntensityRules .
<code>rates</code>	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
<code>intrequired</code>	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
<code>coelCutoff</code>	coelution score threshold between peaks (adducts, parent and fragment ions...). Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>db</code>	list of data bases required for annotation. By default, <code>db</code> contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>db</code> may need to be supplied. See createLipidDB and assignDB .

Details

`idCerpos` function involves 5 steps. 1) FullMS-based identification of candidate Cer as M+H, M+H-H₂O and M+Na. 2) Search of Cer class fragments: there isn't any class fragment by default. 3) Search of specific fragments that inform about the sphingoid base (Sph as M+H-2H₂O resulting from the loss of the FA chain) and the FA chain (by default it is calculated using the difference between precursor and sph fragments). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In this case, there are no intensity rules by default.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (*m/z* error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `msobject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity), Adducts, ppm (*m/z* error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafes.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idCerpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idCLneg

Cardiolipines (CL) annotation for ESI-

Description

CL identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idCLneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 5,
  rt,
  adducts = c("M-H", "M+Na-2H"),
  clfrags = c(),
  clrequired = c(),
  ftype = c(),
  chainfrags_sn1 = c("lysopa_M-H-H2O"),
  chainfrags_sn2 = c("lysopa_M-H-H2O"),
  chainfrags_sn3 = c("lysopa_M-H-H2O"),
  chainfrags_sn4 = c("lysopa_M-H-H2O"),
  intrules = c("Unknown"),
  rates = c(),
  intrequired = c(),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

mobject	an mobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursor and product ions. By default, 3 seconds.
rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for CL in ESI-. Adducts allowed can be modified in adductsTable (dbs argument).
clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
chainfrags_sn1	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFragments for details.
chainfrags_sn2	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFragments for details.
chainfrags_sn3	character vector containing the fragmentation rules for the chain fragments in sn3 position. See chainFragments for details.
chainfrags_sn4	character vector containing the fragmentation rules for the chain fragments in sn4 position. See chainFragments for details.
intrules	character vector specifying the fragments to compare. See checkIntensityRules . If some intensity rules should be employed to identify the chains position but they are't known yet, use "Unknown". If it isn't required, leave an empty vector.
rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idCLneg function involves 5 steps. 1) FullMS-based identification of candidate CL as M-H or M-2H. 2) Search of CL class fragments: no class fragments are searched by defaults as they use to have bad coelution scores. 3) Search of specific fragments that inform about chain composition at

sn1 (lysoPA as M-H-H₂O), sn2 (lysoPA as M-H-H₂O), sn3 (lysoPA as M-H-H₂O) and sn4 (lysoPA as M-H-H₂O). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. For CL there are no intensity rules by default.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idCLneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idDGpos

Diacylglycerols (DG) annotation for ESI+

Description

DG identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idDGpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H-H2O", "M+NH4", "M+Na"),
  clfrags = c(),
  clrequired = c(),
  ftype = c(),
  chainfrags_sn1 = c("mg_M+H-H2O"),
  chainfrags_sn2 = c("mg_M+H-H2O"),
  intrules = c("mg_sn1/mg_sn2"),
  rates = c("1"),
  intrequired = c(T),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>msobject</code>	an <code>msobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for DG in ESI+. Adducts allowed can be modified in <code>adductsTable</code> (<code>dbs</code> argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments in <code>sn1</code> position. See chainFragments for details.
<code>chainfrags_sn2</code>	character vector containing the fragmentation rules for the chain fragments in <code>sn2</code> position. See chainFragments for details. If empty, it will be estimated based on the difference between precursors and <code>sn1</code> chains.
<code>intrules</code>	character vector specifying the fragments to compare. See checkIntensityRules .

rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
db	list of data bases required for annotation. By default, db contains the required data frames based on the default fragmentation rules. If these rules are modified, db may need to be supplied. See createLipidDB and assignDB .

Details

idDGpos function involves 5 steps. 1) FullMS-based identification of candidate DG as M+H-H₂O, M+NH₄ and M+Na. 2) Search of DG class fragments: there are no class fragment by default. 3) Search of specific fragments that inform about the FA chains (MGs as M+H-H₂O resulting from the loss of the FA chains). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position: MG coming from the loss of the sn2 chain is more intense than the one coming from the loss of sn1.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m/z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity), Adducts, ppm (m/z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")
```



```
library(LipidMS)
msobject <- idDGpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idFAHFAneg

FAHFA annotation for ESI-

Description

FAHFA identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idFAHFAneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H"),
  clfrags = c(),
  clrequired = c(),
  ftype = c(),
  chainfrags_sn1 = c("hfa_M-H"),
  chainfrags_sn2 = c("fa_M-H"),
  intrules = c("hfa_sn1/fa_sn2"),
  rates = c("3/1"),
  intrequired = c(T),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

msobject	an msobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursor and product ions. By default, 3 seconds.
rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for FAHFA in ESI-. Adducts allowed can be modified in adductsTable (dbs argument).

clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
chainfrags_sn1	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFragments for details.
chainfrags_sn2	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFragments for details. If empty, it will be estimated based on the difference between precursors and sn1 chains.
intrules	character vector specifying the fragments to compare. See checkIntensityRules .
rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
db	list of data bases required for annotation. By default, db contains the required data frames based on the default fragmentation rules. If these rules are modified, db may need to be supplied. See createLipidDB and assignDB .

Details

idFAHFAneg function involves 5 steps. 1) FullMS-based identification of candidate FAHFA as M-H. 2) Search of FAHFA class fragments: there isn't any class fragment by default. 3) Search of specific fragments that inform about chain composition in sn1 (HFA as M-H resulting from the loss of the FA chain) and sn2 (FA chain as M-H). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In this case, HFA intensity has to be higher than FA.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idFAHFAneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idFAneg

Fatty Acids (FA) annotation for ESI-

Description

FA identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idFAneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H", "2M-H"),
  clfrags = c("fa_M-H", "fa_M-H-H2O"),
  clrequired = c(FALSE, FALSE),
  ftype = c("BB", "BB"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>msoject</code>	an <code>msoject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for FA in ESI-. Adducts allowed can be modified in <code>adductsTable</code> (<code>db</code> s argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>db</code> s	list of data bases required for annotation. By default, <code>db</code> s contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>db</code> s may need to be supplied. See createLipidDB and assignDB .

Details

`idFAneg` function involves 2 steps. 1) FullMS-based identification of candidate FA as M-H or 2M-H. 2) Search of FA class fragments: neutral loss of H₂O coeluting with the precursor ion or the molecular ion.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (*m/z* error), confidenceLevel (in this case, just MS-only or Subclass level (if any class fragment is defined) are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `msoject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity), Adducts, ppm (*m/z* error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the `annotatedPeaklist` element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idFAneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idLPCneg*Lysophosphocholines (LPC) annotation for ESI-*

Description

LPC identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idLPCneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+CH3COO", "M-CH3", "M+CH3COO-CH3"),
  clfrags = c(168.0426, 224.0688, "lysopa_M-H", "lysopc_M-CH3"),
  clrequired = c(F, F, F, F),
  ftype = c("F", "F", "BB", "BB"),
  chainfrags_sn1 = c("fa_M-H"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>mobject</code>	an <code>mobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for LPC in ESI-. Adducts allowed can be modified in <code>adductsTable</code> (<code>db</code> s argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments. See chainFrgs for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>db</code> s	list of data bases required for annotation. By default, <code>db</code> s contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>db</code> s may need to be supplied. See createLipidDB and assignDB .

Details

`idLPCneg` function involves 3 steps. 1) FullMS-based identification of candidate LPC as $M+CH_3COO$, $M-CH_3$ and $M+CH_3COO-CH_3$. To avoid incorrect annotations of PE as PC, candidates which are present just as $M-CH_3$ will be ignored. 2) Search of LPC class fragments: 168.0426, 224.0688, lysoPA as $M-H$ or lysoPC as $M-CH_3$ coeluting with the precursor ion. 3) Search of specific fragments that confirm chain composition (FA as $M-H$).

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (*m/z* error), `confidenceLevel` (in this case, as LPC only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `mobject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity), Adducts, ppm (*m/z* error), `confidenceLevel` (Subclass, FA level, where chains are known but not their positions, or FA position level), `peakID`, and PFCS (parent-fragment coelution score mean of all fragments used

for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idLPCneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idLPCpos

Lysophosphocholines (LPC) annotation for ESI+

Description

LPC identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idLPCpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H", "M+Na"),
  clfrags = c(104.1075, 184.0739),
  clrequired = c(F, F),
  ftype = c("F", "F"),
  chainfrags_sn1 = c("mg_M+H-H2O"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>mobject</code>	an <code>mobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for LPC in ESI+. Adducts allowed can be modified in <code>adductsTable</code> (<code>db</code> s argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments. See chainFragments for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>db</code> s	list of data bases required for annotation. By default, <code>db</code> s contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>db</code> s may need to be supplied. See createLipidDB and assignDB .

Details

`idLPCpos` function involves 3 steps. 1) FullMS-based identification of candidate LPC as $M+H$ and $M+Na$. 2) Search of LPC class fragments: 104.1075 and 184.0739 coeluting with the precursor ion. 3) Search of specific fragments that confirm chain composition (MG as $M+H-H_2O$).

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity, which comes directly from `de` input), Adducts, ppm (*m/z* error), `confidenceLevel` (in this case, as LPC only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `mobject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity), Adducts, ppm (*m/z* error), `confidenceLevel` (Subclass, FA level, where chains are known but not their positions, or FA position level), `peakID`, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the `annotatedPeaklist` element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idLPCpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idLPEneg

Lysophosphoethanolamines (LPE) annotation for ESI-

Description

LPE identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idLPEneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H"),
  clfrags = c(140.0115, 196.038, 214.048, "lysope_M-CH3"),
  clrequired = c(F, F, F, "excluding"),
  ftype = c("F", "F", "F", "BB"),
  chainfrags_sn1 = c("fa_M-H"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>mobject</code>	an <code>mobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for LPE in ESI-. Adducts allowed can be modified in <code>adductsTable</code> (<code>db</code> s argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments. See chainFragments for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>db</code> s	list of data bases required for annotation. By default, <code>db</code> s contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>db</code> s may need to be supplied. See createLipidDB and assignDB .

Details

`idLPEneg` function involves 3 steps. 1) FullMS-based identification of candidate LPE as M-H. 2) Search of LPE class fragments: 140.0115, 196.038 and 214.048 coeluting with the precursor ion. If a loss of CH₃ group is found coeluting with any candidate, this will be excluded as it is a characteristic fragment of LPC. 3) Search of specific fragments that confirm chain composition (FA as M-H).

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m/z error), confidenceLevel (in this case, as LPE only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `mobject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity), Adducts, ppm (m/z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used

for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idLPEneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idLPEpos

Lysophosphoethanolamines (LPE) annotation for ESI+

Description

LPE identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idLPEpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H", "M+Na"),
  clfrags = c(141.01909),
  clrequired = c(F),
  ftype = c("NL"),
  chainfrags_sn1 = c("mg_M+H-H2O"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>mobject</code>	an <code>mobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for LPE in ESI+. Adducts allowed can be modified in <code>adductsTable</code> (<code>db</code> s argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments. See chainFragments for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>db</code> s	list of data bases required for annotation. By default, <code>db</code> s contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>db</code> s may need to be supplied. See createLipidDB and assignDB .

Details

`idLPEpos` function involves 3 steps. 1) FullMS-based identification of candidate LPE as $M+H$ and $M+Na$. 2) Search of LPE class fragments: neutral loss of 141.01909 coeluting with the precursor ion. 3) Search of specific fragments that confirm chain composition in `sn1` (MG as $M+H-H_2O$).

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), `mz`, RT (in seconds), I (intensity, which comes directly from `de` input), Adducts, ppm (`m.z` error), `confidenceLevel` (in this case, as LPE only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `mobject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), `m.z`, RT (in seconds), I (intensity), Adducts, ppm (`m.z` error), `confidenceLevel` (Subclass, FA level, where chains are known but not their positions, or FA position level), `peakID`, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the `annotatedPeaklist` element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idLPEpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idLPGneg

Lysophosphoglycerols (LPG) annotation for ESI-

Description

LPG identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idLPGneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H"),
  clfrags = c(152.9958, 227.0326, 209.022, 74.0359),
  clrequired = c(F, F, F, F),
  ftype = c("F", "F", "F", "NL"),
  chainfrags_sn1 = c("fa_M-H"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>mobject</code>	an <code>mobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for LPG in ESI-. Adducts allowed can be modified in <code>adductsTable</code> (<code>db</code> s argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments. See chainFragments for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>db</code> s	list of data bases required for annotation. By default, <code>db</code> s contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>db</code> s may need to be supplied. See createLipidDB and assignDB .

Details

`idLPGneg` function involves 3 steps. 1) FullMS-based identification of candidate LPG as M-H. 2) Search of LPG class fragments: 152.9958, 227.0326, 209.022 and neutral loss of 74.0359 coeluting with the precursor ion. 3) Search of specific fragments that confirm chain composition (FA as M-H).

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity, which comes directly from `de` input), Adducts, ppm (*m/z* error), `confidenceLevel` (in this case, as LPG only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `mobject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity), Adducts, ppm (*m/z* error), `confidenceLevel` (Subclass, FA level, where chains are known but not their positions, or FA position level), `peakID`, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the `annotatedPeaklist` element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idLPGneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idLPIneg

Lysophosphoinositols (LPI) annotation for ESI-

Description

LPI identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idLPIneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H"),
  clfrags = c(241.0115, 223.0008, 259.0219, 297.0375),
  clrequired = c(F, F, F, F),
  ftype = c("F", "F", "F", "F"),
  chainfrags_sn1 = c("fa_M-H"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>mobject</code>	an <code>mobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for LPI in ESI-. Adducts allowed can be modified in <code>adductsTable</code> (<code>db</code> s argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments. See chainFragments for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>db</code> s	list of data bases required for annotation. By default, <code>db</code> s contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>db</code> s may need to be supplied. See createLipidDB and assignDB .

Details

`idLPIneg` function involves 3 steps. 1) FullMS-based identification of candidate LPI as M-H. 2) Search of LPI class fragments: 241.0115, 223.0008, 259.0219 and 297.0375 coeluting with the precursor ion. 3) Search of specific fragments that confirm chain composition (FA as M-H).

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity, which comes directly from `de` input), Adducts, ppm (*m/z* error), `confidenceLevel` (in this case, as LPI only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `mobject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity), Adducts, ppm (*m/z* error), `confidenceLevel` (Subclass, FA level, where chains are known but not their positions, or FA position level), `peakID`, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the `annotatedPeaklist` element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idLPIneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idLPSneg

Lysophosphoserines (LPS) annotation for ESI-

Description

LPS identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idLPSneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H", "M+Na-2H"),
  clfrags = c(87.032),
  clrequired = c(F),
  ftype = c("NL"),
  chainfrags_sn1 = c("fa_M-H"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>mobject</code>	an <code>mobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for LPS in ESI-. Adducts allowed can be modified in <code>adductsTable</code> (<code>db</code> s argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments. See chainFragments for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>db</code> s	list of data bases required for annotation. By default, <code>db</code> s contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>db</code> s may need to be supplied. See createLipidDB and assignDB .

Details

`idLPSneg` function involves 3 steps. 1) FullMS-based identification of candidate LPS as M-H and M+Na-2H. 2) Search of LPS class fragments: neutral loss of 87.032 coeluting with the precursor ion. 3) Search of specific fragments that confirm chain composition (FA as M-H).

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m/z error), confidenceLevel (in this case, as LPS only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `mobject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity), Adducts, ppm (m/z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the `annotatedPeaklist` element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idLPSneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idMGpos

Monoacylglycerol (MG) annotation for ESI+

Description

MG identification based on fragmentation patterns for LC-MS/MS DIA and DDA data acquired in positive mode.

Usage

```
idMGpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H-H2O", "M+NH4", "M+Na"),
  clfrags = c(),
  clrequired = c(),
  ftype = c(),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>mobject</code>	an <code>mobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for MG in ESI+. Adducts allowed can be modified in <code>adductsTable</code> (<code>db</code> s argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>db</code> s	list of data bases required for annotation. By default, <code>db</code> s contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>db</code> s may need to be supplied. See createLipidDB and assignDB .

Details

`idMGpos` function involves 2 steps. 1) FullMS-based identification of candidate MG as $M+H-H_2O$, $M+NH_4$ and $M+Na$. 2) Search of MG class fragments if any is assigned.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (*m/z* error), confidenceLevel (in this case, just MS-only or Subclass level (if any class fragment is defined) are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `mobject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), *m/z*, RT (in seconds), I (intensity), Adducts, ppm (*m/z* error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the `annotatedPeaklist` element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idMGpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idNEG

Lipids annotation for ESI-

Description

Lipids annotation based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode. This function compiles all functions written for ESI- annotations.

Usage

```
idNEG(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 5,
  coelCutoff = 0.8,
  lipidClasses = c("FA", "FAHFA", "LPC", "LPE", "LPG", "LPI", "LPS", "PC", "PE", "PG",
    "PI", "PS", "Sph", "SphP", "Cer", "CL", "BA"),
  dbs
)
```

Arguments

msobject	an msobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.

rttol	total rt window for coelution between precursor and product ions. By default, 5 seconds.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
lipidClasses	classes of interest to run the identification functions.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idNEG(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idPCneg

Phosphocholines (PC) annotation for ESI-

Description

PC identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idPCneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+CH3COO", "M-CH3", "M+CH3COO-CH3"),
  clfrags = c(168.0426, 224.0688, "pc_M-CH3"),
  clrequired = c(F, F, F),
  ftype = c("F", "F", "BB"),
  chainfrags_sn1 = c("lysopc_M-CH3"),
  chainfrags_sn2 = c("fa_M-H", "lysopc_M-CH3"),
  intrules = c("lysopc_sn1/lysopc_sn2"),
  rates = c("3/1"),
  intrequired = c(T),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>msobject</code>	an <code>msobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for PC in ESI-. Adducts allowed can be modified in <code>adductsTable</code> (<code>dbs</code> argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments in <code>sn1</code> position. See chainFragments for details.
<code>chainfrags_sn2</code>	character vector containing the fragmentation rules for the chain fragments in <code>sn2</code> position. See chainFragments for details. If empty, it will be estimated based on the difference between precursors and <code>sn1</code> chains.
<code>intrules</code>	character vector specifying the fragments to compare. See checkIntensityRules .

rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
db	list of data bases required for annotation. By default, db contains the required data frames based on the default fragmentation rules. If these rules are modified, db may need to be supplied. See createLipidDB and assignDB .

Details

idPCneg function involves 5 steps. 1) FullMS-based identification of candidate PC as M+CH₃COO, M-CH₃ or M+CH₃COO-CH₃. To avoid incorrect annotations of PE as PC, candidates which are present just as M-CH₃ will be ignored. 2) Search of PC class fragments: 168.0426, 224.0688 or loss of CH₃ coeluting with the precursor ion. 3) Search of specific fragments that inform about chain composition in sn1 (lysoPC as M-CH₃ resulting from the loss of the FA chain at sn2) and sn2 (lysoPC as M-CH₃ resulting from the loss of sn1 or FA as M-H). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In this case, lysoPC from sn1 is at least 3 times more intense than lysoPC from sn2.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m/z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity), Adducts, ppm (m/z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPCneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idPCpos	<i>Phosphocholines (PC) annotation for ESI+</i>
---------	---

Description

PC identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idPCpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H", "M+Na"),
  clfrags = c(104.1075, 184.0739, 183.06604),
  clrequired = c(F, F, F),
  ftype = c("F", "F", "NL"),
  chainfrags_sn1 = c("lysopc_M+H", "lysopc_M+H-H2O"),
  chainfrags_sn2 = c("lysopc_M+H", "lysopc_M+H-H2O", ""),
  intrules = c("lysopc_sn1/lysopc_sn2"),
  rates = c("2/1"),
  intrequired = c(T),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

msobject	an msobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursor and product ions. By default, 3 seconds.

rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for PC in ESI+. Adducts allowed can be modified in adductsTable (dbs argument).
clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
chainfrags_sn1	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFrgs for details.
chainfrags_sn2	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFrgs for details. If empty, it will be estimated based on the difference between precursors and sn1 chains.
intrules	character vector specifying the fragments to compare. See checkIntensityRules .
rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idPCpos function involves 5 steps. 1) FullMS-based identification of candidate PC as M+H and M+Na. 2) Search of PC class fragments: 104.1075, 184.0739 and neutral loss of 183.06604 coeluting with the precursor ion. 3) Search of specific fragments that inform about chain composition in sn1 (lysoPC as M+H or M+H-H₂O resulting from the loss of the FA chain at sn2) and sn2 (lysoPC as M+H or M+H-H₂O resulting from the loss of the FA chain at sn1 or the difference between precursor and sn1 chain fragments). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In this case, lysoPC from sn1 is at least twice more intense than lysoPC from sn2.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPCpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idPEneg

Phosphoethanolamines (PE) annotation for ESI-

Description

PE identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idPEneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 5,
  rt,
  adducts = c("M-H"),
  clfrags = c(140.0118, 196.038, 214.048, "pe_M-CH3"),
```

```

    clrequired = c(F, F, F, "excluding"),
    ftype = c("F", "F", "F", "BB"),
    chainfrags_sn1 = c("lysope_M-H"),
    chainfrags_sn2 = c("lysope_M-H", "fa_M-H"),
    intrules = c("lysope_sn1/lysope_sn2"),
    rates = c("3/1"),
    intrequired = c(T),
    coelCutoff = 0.8,
    dbs
  )

```

Arguments

<code>msobject</code>	an <code>msobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for PE in ESI-. Adducts allowed can be modified in <code>adductsTable</code> (<code>dbs</code> argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFragments for details.
<code>chainfrags_sn2</code>	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFragments for details. If empty, it will be estimated based on the difference between precursors and sn1 chains.
<code>intrules</code>	character vector specifying the fragments to compare. See checkIntensityRules .
<code>rates</code>	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
<code>intrequired</code>	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>dbs</code>	list of data bases required for annotation. By default, <code>dbs</code> contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>dbs</code> may need to be supplied. See createLipidDB and assignDB .

Details

idPEneg function involves 5 steps. 1) FullMS-based identification of candidate PE as M-H. 2) Search of PE class fragments: 140.0115, 196.038, 214.048 ion coeluting with the precursor ion. If a loss of CH₃ group is found coeluting with any candidate, this will be excluded as it is a characteristic fragment of PC. 3) Search of specific fragments that inform about chain composition in sn1 (lysoPE as M-H resulting from the loss of the FA chain at sn2) and sn2 (lysoPE as M-H resulting from the loss of the FA chain at sn1 or FA chain as M-H). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In this case, lysoPE from sn1 is at least 3 times more intense than lysoPE from sn2.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPEneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idPEpos

*Phosphoethanolamines (PE) annotation for ESI+***Description**

PE identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idPEpos(
  msubject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H", "M+Na"),
  clfrags = c("dg_M+H-H2O"),
  clrequired = c(F),
  ftype = c("BB"),
  chainfrags_sn1 = c("lysope_M+H-H2O", "mg_M+H-H2O"),
  chainfrags_sn2 = c("mg_M+H-H2O"),
  intrules = c("lysope_sn1/lysope_sn1", "mg_sn1/mg_sn2"),
  rates = c("3/1", "1/2"),
  intrequired = c(F, F),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>msubject</code>	an <code>msobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total <code>rt</code> window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	<code>rt</code> range where the function will look for candidates. By default, it will search within all <code>RT</code> range in <code>MS1</code> .
<code>adducts</code>	expected adducts for PE in ESI+. Adducts allowed can be modified in <code>adductsTable</code> (<code>dbs</code> argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.

ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
chainfrags_sn1	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFrag s for details.
chainfrags_sn2	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFrag s for details. If empty, it will be estimated based on the difference between precursors and sn1 chains.
intrules	character vector specifying the fragments to compare. See checkIntensityRules .
rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idPEpos function involves 5 steps. 1) FullMS-based identification of candidate PE as M+H and M+Na. 2) Search of PE class fragments: loss of head group (DG as M+H-H₂O) coeluting with the precursor ion. 3) Search of specific fragments that inform about chain composition at sn1 (MG as M+H-H₂O resulting from the loss of the FA chain at sn2 and the head group or LPE as M+H-H₂O resulting just from the loss of the FA chain) and sn2 (FA or MG chain from sn2 as M+H-H₂O or the difference between precursor and sn1 chain fragments). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. LPE or MG from sn1 is at least 3 times more intense than the ones from sn2.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafes.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPEpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idPGneg

Phosphoglycerols (PG) annotation for ESI-

Description

PG identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idPGneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H"),
  clfrags = c(152.9958, 227.0326, 209.022, 74.0359),
  clrequired = c(F, F, F, F),
  ftype = c("F", "F", "F", "NL"),
  chainfrags_sn1 = c("lysopg_M-H"),
  chainfrags_sn2 = c("lysopg_M-H", "fa_M-H"),
  intrules = c("lysopg_sn1/lysopg_sn2"),
  rates = c("2/1"),
  intrequired = c(T),
  coelCutoff = 0.8,
  dbs
)
```


Arguments

mobject	an mobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursor and product ions. By default, 3 seconds.
rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for PG in ESI-. Adducts allowed can be modified in adductsTable (dbs argument).
clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
chainfrags_sn1	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFragments for details.
chainfrags_sn2	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFragments for details. If empty, it will be estimated based on the difference between precursors and sn1 chains.
intrules	character vector specifying the fragments to compare. See checkIntensityRules .
rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idPGneg function involves 5 steps. 1) FullMS-based identification of candidate PG as M-H. 2) Search of PG class fragments: 152.9958, 227.0326, 209.022 and neutral loss of 74.0359 coeluting with the precursor ion. 3) Search of specific fragments that inform about chain composition at sn1 (lysoPG as M-H resulting from the loss of the FA chain at sn2) and sn2 (lysoPG as M-H resulting from the loss of the FA chain at sn1 or FA chain as M-H). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In this case, lysoPG from sn1 is at least 3 times more intense than lysoPG from sn2.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity,

which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPGneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idPGpos

Phosphoglycerols (PG) annotation for ESI+

Description

PG identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idPGpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H", "M+NH4", "M+Na"),
  clfrags = c("dg_M+H-H2O"),
  clrequired = c(F),
  ftype = c("BB"),
  chainfrags_sn1 = c("mg_M+H-H2O"),
  chainfrags_sn2 = c("mg_M+H-H2O"),
  intrules = c("mg_sn1/mg_sn2"),
  rates = c("1/2"),
  intrequired = c(F, F),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

<code>msobject</code>	an <code>msobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for PE in ESI+. Adducts allowed can be modified in <code>adductsTable</code> (<code>dbs</code> argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments in <code>sn1</code> position. See chainFragments for details.
<code>chainfrags_sn2</code>	character vector containing the fragmentation rules for the chain fragments in <code>sn2</code> position. See chainFragments for details. If empty, it will be estimated based on the difference between precursors and <code>sn1</code> chains.
<code>intrules</code>	character vector specifying the fragments to compare. See checkIntensityRules .

rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idPGpos function involves 5 steps. 1) FullMS-based identification of candidate PG as M+H, M+NH₄ and M+Na. 2) Search of PG class fragments: loss of head group (DG as M+H-H₂O) coeluting with the precursor ion. 3) Search of specific fragments that inform about chain composition at sn1 (MG as M+H-H₂O resulting from the loss of the FA chain at sn2) and sn2 (MG as M+H-H₂O resulting from the loss of the FA chain at sn1). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. MG from sn1 is at least 2 times more intense than the ones from sn2.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPGpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idPIneg

Phosphoinositols (PI) annotation for ESI-

Description

PI identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idPIneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H"),
  clfrags = c(241.0115, 223.0008, 259.0219, 297.0375),
  clrequired = c(F, F, F, F),
  ftype = c("F", "F", "F", "F"),
  chainfrags_sn1 = c("lysopi_M-H", "lysopa_M-H"),
  chainfrags_sn2 = c("lysopi_M-H", "lysopa_M-H", "fa_M-H"),
  intrules = c("lysopi_sn1/lysopi_sn2", "lysopa_sn1/lysopa_sn2"),
  rates = c("3/1", "3/1"),
  intrequired = c(F, F),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

msobject	an msobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursor and product ions. By default, 3 seconds.

rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for PI in ESI-. Adducts allowed can be modified in adductsTable (dbs argument).
clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
chainfrags_sn1	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFrgs for details.
chainfrags_sn2	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFrgs for details. If empty, it will be estimated based on the difference between precursors and sn1 chains.
intrules	character vector specifying the fragments to compare. See checkIntensityRules .
rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idPIneg function involves 5 steps. 1) FullMS-based identification of candidate PI as M-H. 2) Search of PI class fragments: 241.0115, 223.0008, 259.0219 and 297.0375 coeluting with the precursor ion. 3) Search of specific fragments that inform about chain composition at sn1 (lysoPI as M-H resulting from the loss of the FA chain at sn2 or lysoPA as M-H if it also losses the head group) and sn2 (lysoPI or lysoPA as M-H resulting from the loss of the FA chain at sn1 or FA chain as M-H). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In this case, lysoPI or lysoPA from sn1 is at least 3 times more intense than lysoPI or lysoPA from sn2.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPIneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idPOS

Lipids annotation for ESI+

Description

Lipids annotation based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode. This function compiles all functions written for ESI+ annotations.

Usage

```
idPOS(
  msubject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 5,
  coelCutoff = 0.8,
  lipidClasses = c("MG", "LPC", "LPE", "PC", "PE", "PG", "Sph", "SphP", "Cer", "SM",
    "Car", "CE", "DG", "TG"),
```

```

    dbs
  )

```

Arguments

<code>mobject</code>	an mobject returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 5 seconds.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>lipidClasses</code>	classes of interest to run the identification functions.
<code>dbs</code>	list of data bases required for annotation. By default, <code>dbs</code> contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>dbs</code> may need to be supplied. See createLipidDB and assignDB .

Value

annotated mobject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```

## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
mobject <- idPOS(LipidMSdata2::mobjectDIApos)

## End(Not run)

```


idPSneg

*Phosphoserines (PS) annotation for ESI-***Description**

PS identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idPSneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H", "M+Na-2H"),
  clfrags = c(87.032, 152.9958),
  clrequired = c(F, F),
  ftype = c("NL", "F"),
  chainfrags_sn1 = c("lysopa_M-H", "lysopa_M-H-H2O"),
  chainfrags_sn2 = c("lysopa_M-H", "lysopa_M-H-H2O", "fa_M-H"),
  intrules = c("lysopa_sn1/lysopa_sn2"),
  rates = c("3/1"),
  intrequired = c(T),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

msobject	an msobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursor and product ions. By default, 3 seconds.
rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for PS in ESI-. Adducts allowed can be modified in adductsTable (dbs argument).
clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.

<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFrgs for details.
<code>chainfrags_sn2</code>	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFrgs for details. If empty, it will be estimated based on the difference between precursors and sn1 chains.
<code>intrules</code>	character vector specifying the fragments to compare. See checkIntensityRules .
<code>rates</code>	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
<code>intrequired</code>	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
<code>coelCutoff</code>	coelution score threshold between parent and fragment ions. Only applied if <code>rawData</code> info is supplied. By default, 0.8.
<code>dbs</code>	list of data bases required for annotation. By default, <code>dbs</code> contains the required data frames based on the default fragmentation rules. If these rules are modified, <code>dbs</code> may need to be supplied. See createLipidDB and assignDB .

Details

`idPSneg` function involves 5 steps. 1) FullMS-based identification of candidate PS as M-H or M+Na-2H. 2) Search of PS class fragments: neutral loss of 87.032 (serine) coeluting with the precursor ion. 3) Search of specific fragments that inform about chain composition at sn1 (lysoPA as M-H or M-H-H₂O resulting from the loss of the FA chain at sn2 and the head group) and sn2 (lysoPA as M-H or M-H-H₂O resulting from the loss of the FA chain at sn1 and the head group or FA chain as M-H). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In this case, lysoPA from sn1 is at least 3 times more intense than lysoPA from sn2.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m/z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated `msoject` (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m/z, RT (in seconds), I (intensity), Adducts, ppm (m/z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the `annotatedPeaklist` element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPSneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idSMpos

Sphingomyelins (SM) annotation for ESI+

Description

SM identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idSMpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H", "M+Na"),
  clfrags = c(104.1075, 184.0739, 183.06604),
  clrequired = c(F, F, F),
  ftype = c("F", "F", "NL"),
  chainfrags_sn1 = c("sph_M+H-2H2O"),
  chainfrags_sn2 = c(""),
  intrules = c(),
  rates = c(),
  intrequired = c(),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

mobject	an mobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursor and product ions. By default, 3 seconds.
rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for SM in ESI+. Adducts allowed can be modified in adductsTable (dbs argument).
clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
chainfrags_sn1	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFragments for details.
chainfrags_sn2	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFragments for details. If empty, it will be estimated based on the difference between precursors and sn1 chains.
intrules	character vector specifying the fragments to compare. See checkIntensityRules .
rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idSMpos function involves 5 steps. 1) FullMS-based identification of candidate SM as M+H and M+Na. 2) Search of SM class fragments: 104.1075, 184.0739 and neutral loss of 183.06604 coeluting with the precursor ion. 3) Search of specific fragments that inform about the composition of the sphingoid base (Sph as M+H-2H₂O resulting from the loss of the FA chain) and the FA chain (by default it is calculated using the difference between precursor and sph chain fragments). 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In this case, there are no intensity rules by default as FA chain is unlikely to be detected.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idSMpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idSphneg

Sphingoid bases (Sph) annotation for ESI-

Description

Sph identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idSphneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H"),
  clfrags = c("sph_M-H-H2O", "sph_M-H-2H2O"),
  clrequired = c(F, F),
  ftype = c("BB", "BB"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

msobject	an msobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursor and product ions. By default, 3 seconds.
rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for Sph in ESI-. Adducts allowed can be modified in adductsTable (dbs argument).
clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idSphneg function involves 2 steps. 1) FullMS-based identification of candidate Sph as M-H. 2) Search of Sph class fragments: neutral loss of 1 or 2 H₂O molecules.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity,

which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (in this case, as Sph only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idSphneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idSphPneg

Sphingoid bases phosphate (SphP) annotation for ESI-

Description

SphP identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in negative mode.

Usage

```
idSphPneg(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M-H"),
  clfrags = c(78.9585, 96.9691, "sphP_M-H-H2O"),
  clrequired = c(F, F, F),
  ftype = c("F", "F", "BB"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

msobject	an msobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursor and product ions. By default, 3 seconds.
rt	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for SphP in ESI-. Adducts allowed can be modified in adductsTable (dbs argument).
clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idSphpos function involves 2 steps. 1) FullMS-based identification of candidate SphP as M-H. 2) Search of SphP class fragments: 78.9585, 96.969 or neutral loss of 1 H₂O molecule.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity,

which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (in this case, as SphP only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maialba@alumni.uv.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idSphPneg(LipidMSdata2::msobjectDIAneg)

## End(Not run)
```

idSphpos

Sphingoid bases (Sph) annotation for ESI-

Description

Sph identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idSphpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H"),
  clfrags = c("sph_M+H-H2O", "sph_M+H-2H2O"),
  clrequired = c(F, F),
  ftype = c("BB", "BB"),
  coelCutoff = 0.8,
  dbs
)
```

Arguments

msobject	an msobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursors and product ions. By default, 3 seconds.
rt	rt window where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for Sph in ESI+. Adducts allowed can be modified in adductsTable (dbs argument).
clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idSphpos function involves 2 steps. 1) FullMS-based identification of candidate Sph as M+H. 2) Search of Sph class fragments: neutral loss of 1 or 2 H₂O molecules.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity,

which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (in this case, as Sph only have one chain, only Subclass and FA level are possible) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idSphpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idSphPpos

Sphingoid bases phosphate (SphP) annotation for ESI+

Description

SphP identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```
idSphPpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+H"),
  clfrags = c("sphP_M+H-H2O", "sphP_M+H-2H2O", "sphP_M+H-H2O-NH4"),
  clrequired = c(F, F, F),
  ftype = c("BB", "BB", "BB"),
  coelCutoff = 0.7,
  dbs
)
```

Arguments

msobject	an msobject returned by dataProcessing .
ppm_precursor	mass tolerance for precursor ions. By default, 5 ppm.
ppm_products	mass tolerance for product ions. By default, 10 ppm.
rttol	total rt window for coelution between precursors and product ions. By default, 3 seconds.
rt	rt window where the function will look for candidates. By default, it will search within all RT range in MS1.
adducts	expected adducts for Sph in ESI+. Adducts allowed can be modified in adductsTable (dbs argument).
clfrags	vector containing the expected fragments for a given lipid class. See checkClass for details.
clrequired	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
ftype	character vector indicating the type of fragments in clfrags. It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idSphPpos function involves 2 steps. 1) FullMS-based identification of candidate SphP as M+H. 2) Search of SphP class fragments: neutral loss of 1 or 2 H₂O molecules, or H₂O and NH₄.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), mz, RT (in seconds), I (intensity,

which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (in this case, as SphP only have one chain, only Subclass and FA level are possible). and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idSphPpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

idTGpos

Triacylglycerols (TG) annotation for ESI+

Description

TG identification based on fragmentation patterns for LC-MS/MS DIA or DDA data acquired in positive mode.

Usage

```

idTGpos(
  msobject,
  ppm_precursor = 5,
  ppm_products = 10,
  rttol = 3,
  rt,
  adducts = c("M+NH4", "M+Na"),
  clfrags = c(),
  clrequired = c(),
  ftype = c(),
  chainfrags_sn1 = c("cbdiff-dg_M+H-H20"),
  chainfrags_sn2 = c("cbdiff-dg_M+H-H20"),
  chainfrags_sn3 = c("cbdiff-dg_M+H-H20"),
  intrules = c("cbdiff-dg_sn2/cbdiff-dg_sn1", "cbdiff-dg_sn2/cbdiff-dg_sn3",
    "cbdiff-dg_sn1/cbdiff-dg_sn3"),
  rates = c("1", "1", "1"),
  intrequired = c(T, T, T),
  coelCutoff = 0.8,
  dbs
)

```

Arguments

<code>msobject</code>	an <code>msobject</code> returned by dataProcessing .
<code>ppm_precursor</code>	mass tolerance for precursor ions. By default, 5 ppm.
<code>ppm_products</code>	mass tolerance for product ions. By default, 10 ppm.
<code>rttol</code>	total rt window for coelution between precursor and product ions. By default, 3 seconds.
<code>rt</code>	rt range where the function will look for candidates. By default, it will search within all RT range in MS1.
<code>adducts</code>	expected adducts for TG in ESI+. Adducts allowed can be modified in <code>adductsTable</code> (<code>dbs</code> argument).
<code>clfrags</code>	vector containing the expected fragments for a given lipid class. See checkClass for details.
<code>clrequired</code>	logical vector indicating if each class fragment is required or not. If any of them is required, at least one of them must be present within the coeluting fragments. See checkClass for details.
<code>ftype</code>	character vector indicating the type of fragments in <code>clfrags</code> . It can be: "F" (fragment), "NL" (neutral loss) or "BB" (building block). See checkClass for details.
<code>chainfrags_sn1</code>	character vector containing the fragmentation rules for the chain fragments in sn1 position. See chainFragments for details.
<code>chainfrags_sn2</code>	character vector containing the fragmentation rules for the chain fragments in sn2 position. See chainFragments for details. If empty, it will be estimated based on the difference between precursors and sn1 chains.

chainfrags_sn3	character vector containing the fragmentation rules for the chain fragments in sn3 position. See chainFragments for details. If empty, it will be estimated based on the difference between precursors and sn2 chains.
intrules	character vector specifying the fragments to compare. See checkIntensityRules . If some intensity rules should be employed to identify the chains position but they are't known yet, use "Unknown". If it isn't required, leave an empty vector.
rates	character vector with the expected rates between fragments given as a string (e.g. "3/1"). See checkIntensityRules .
intrequired	logical vector indicating if any of the rules is required. If not, at least one must be verified to confirm the structure.
coelCutoff	coelution score threshold between parent and fragment ions. Only applied if rawData info is supplied. By default, 0.8.
dbs	list of data bases required for annotation. By default, dbs contains the required data frames based on the default fragmentation rules. If these rules are modified, dbs may need to be supplied. See createLipidDB and assignDB .

Details

idTGpos function involves 5 steps. 1) FullMS-based identification of candidate TG as M+NH₄ and M+Na. 2) Search of TG class fragments: there are no class fragment by default. 3) Search of specific fragments that inform about the FA chains: DGs resulting from the loss of FA chains as M+H-H₂O. 4) Look for possible chains structure based on the combination of chain fragments. 5) Check intensity rules to confirm chains position. In the case of TG, DG resulting from the loss of sn2 if the most intense, followed by the loss of sn1 and sn3, but this FA position level still needs to be improved due to the high level of coelution for TG.

Results data frame shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity, which comes directly from de input), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level) and PFCS (parent-fragment coelution score mean of all fragments used for the identification).

Value

annotated msubject (list with several elements). The results element is a data frame that shows: ID, class of lipid, CDB (total number of carbons and double bounds), FA composition (specific chains composition if it has been confirmed), m.z, RT (in seconds), I (intensity), Adducts, ppm (m.z error), confidenceLevel (Subclass, FA level, where chains are known but not their positions, or FA position level), peakID, and PFCS (parent-fragment coelution score mean of all fragments used for the identification); and the annotatedPeaklist element shows the original MS1 peaklist with the annotations on it.

Note

This function has been written based on fragmentation patterns observed for three different platforms (QTOF 6550 from Agilent, Sinapt G2-Si from Waters and Q-exactive from Thermo), but it may need to be customized for other platforms or acquisition settings.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idTGpos(LipidMSdata2::msobjectDIApos)

## End(Not run)
```

LipidMSapp

LipidMS shiny app

Description

Interactive UI for LipidMS

Usage

```
LipidMSapp()
```

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:
# example data files can be download from github.com/maialba3/LipidMSv2.0_exampleFiles

library(LipidMS)
LipidMSapp()

## End(Not run)
```

lysopadb

LPA's database

Description

In silico generated database for common LPAs.

Usage

```
data("lysopadb")
```

Format

Data frame with 30 observations and the following 3 variables.

`formula` character vector containing molecular formulas.

`total` character vector indicating the total number of carbons and double bounds of the chains.

`Mass` numeric vector with the neutral masses.

lysopcdb

LPC's database

Description

In silico generated database for common LPCs.

Usage

```
data("lysopcdb")
```

Format

Data frame with 30 observations and the following 3 variables.

`formula` character vector containing molecular formulas.

`total` character vector indicating the total number of carbons and double bounds of the chains.

`Mass` numeric vector with the neutral masses.

lysopedb

LPEs database

Description

In silico generated database for common LPEs.

Usage

```
data("lysopedb")
```

Format

Data frame with 30 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

lysopgdb

LPGs database

Description

In silico generated database for common LPGs.

Usage

```
data("lysopgdb")
```

Format

Data frame with 30 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

lysopidb

LPIs database

Description

In silico generated database for common LPIs.

Usage

```
data("lysopidb")
```

Format

Data frame with 30 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

lysopsdb

LPSs database

Description

In silico generated database for common LPSs

Usage

```
data("lysopsdb")
```

Format

Data frame with 30 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

mgdb

MGs database

Description

In silico generated database for common MGs.

Usage

```
data("mgdb")
```

Format

Data frame with 30 observations and the following 3 variables.

`formula` character vector containing molecular formulas.

`total` character vector indicating the total number of carbons and double bounds of the chains.

`Mass` numeric vector with the neutral masses.

nlsphdb

Neutral losses db for sphingoid bases. It is employed by idCerneg function.

Description

In silico generated database for neutral losses of sphingoid bases in ESI-.

Usage

```
data("nlsphdb")
```

Format

Data frame with 4 observations and the following 3 variables.

`formula` character vector containing molecular formulas.

`total` character vector indicating the total number of carbons and double bounds of the chains.

`Mass` numeric vector with the neutral masses.

organizeResults	<i>Prepare output for LipidMS annotation functions</i>
-----------------	--

Description

Prepare a readable output for LipidMS identification functions.

Usage

```
organizeResults(  
  candidates,  
  clfrags,  
  classConf,  
  chainsComb,  
  intrules,  
  intConf,  
  nchains,  
  class  
)
```

Arguments

candidates	candidates data frame. Output of findCandidates .
clfrags	vector containing the expected fragments for a given lipid class.
classConf	output of checkClass
chainsComb	output of combineChains
intrules	character vector specifying the fragments to compare. See checkIntensityRules .
intConf	output of checkIntensityRules
nchains	number of chains of the targeted lipid class.
class	character value. Lipid class (i.e. PC, PE, DG, TG, etc.).

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

Examples

```
## Not run:  
devtools::install_github("maialba3/LipidMSdata2")  
  
library(LipidMS)  
dbs <- assignDB()  
  
MS1 <- LipidMSdata2::msobjectDIAneg$peaklist$MS1  
MS1 <- MS1[MS1$isotope %in% c("[M+0]"), !colnames(MS1) %in% c("isotope", "group")]  
MS2 <- LipidMSdata2::msobjectDIAneg$peaklist$MS2[,c("m.z", "RT", "int", "peakID")]
```

```

rawData <- rbind(LipidMSdata2::msobjectDIAneg$MS1, LipidMSdata2::msobjectDIAneg$MS2)

candidates <- findCandidates(MS1 = MS1, db = dbs$pgdb, ppm = 10,
rt = c(0, 2000), adducts = c("M-H"), rttol = 10, dbs = dbs,
rawData = rawData, coelCutoff = 0.8)
coelfrags <- coelutingFragments(candidates, MSMS, rttol = 10, rawData = rawData,
coelCutoff = 0.8)
classConf <- checkClass(candidates, coelfrags,
clfrags = c(227.0326, 209.022, 74.0359),
clrequisites = c(FALSE, FALSE, FALSE, FALSE),
ftype = c("F", "F", "NL"), ppm = 10, dbs = dbs)
sn1 <- chainFragments(coelfrags, chainfrags = c("lysopg_M-H"),
candidates = candidates, ppm = 10, dbs = dbs)
sn2 <- chainFragments(coelfrags, chainfrags = c("fa_M-H"), ppm = 10,
candidates = candidates, dbs = dbs)
chainsComb <- combineChains(candidates, nchains=2, sn1, sn2)
intConf <- checkIntensityRules(inrules = c("lysopg_sn1/lysopg_sn1"),
rates = c("2/1"), intrequired = c(TRUE), nchains=2, chainsComb)

res <- organizeResults(candidates, clfrags = c(227.0326, 209.022, 74.0359),
classConf, chainsComb, intrules = c("lysopg_sn1/lysopg_sn1"), intConf,
nchains = 2, class="PG")

## End(Not run)

```

padb

PAs database

Description

In silico generated database for common PAs.

Usage

```
data("padb")
```

Format

Data frame with 147 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

pcdb*PCs database*

Description

In silico generated database for common PCs.

Usage

```
data("pcdb")
```

Format

Data frame with 147 observations and the following 3 variables.

`formula` character vector containing molecular formulas.

`total` character vector indicating the total number of carbons and double bounds of the chains.

`Mass` numeric vector with the neutral masses.

pedb*PEs database*

Description

In silico generated database for common PEs.

Usage

```
data("pedb")
```

Format

Data frame with 147 observations and the following 3 variables.

`formula` character vector containing molecular formulas.

`total` character vector indicating the total number of carbons and double bounds of the chains.

`Mass` numeric vector with the neutral masses.

pgdb

PGs database

Description

In silico generated database for common PGs.

Usage

```
data("pgdb")
```

Format

Data frame with 147 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

pidb

PIs database

Description

In silico generated database for common PIs.

Usage

```
data("pidb")
```

Format

Data frame with 147 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

`plotLipids`*Plot informative peaks for lipids annotation*

Description

Plot informative peaks for each lipid annotated using idPOS and idNEG (or similar functions).

Usage

```
plotLipids(msobject, spar = 0.4)
```

Arguments

<code>msobject</code>	annotated msobject.
<code>spar</code>	smoothing parameter. Numeric value between 0 and 1.

Details

Peak intensities are relative to the maximum intensity of each peak to ease visualization.

Value

`msobject` with a `plots` element which contains a list of plots. Plots on the left side represent raw values while plots on the left side are smoothed or clean scans (MS2 in DDA).

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafes.es>

Examples

```
## Not run:
devtools::install_github("maialba3/LipidMSdata2")

library(LipidMS)
msobject <- idPOS(LipidMSdata2::msobjectDIApos)
msobject <- plotLipids(msobject)

# display the first plot
msobject$plots[[1]]
msobject$plots[["peakID"]]

# save all plot to a pdf file
pdf("plotresults.pdf")
msobject$plots
dev.off()

## End(Not run)
```

psdb

*PSs database***Description**

In silico generated database for common PSs.

Usage

```
data("psdb")
```

Format

Data frame with 147 observations and the following 3 variables.

`formula` character vector containing molecular formulas.

`total` character vector indicating the total number of carbons and double bounds of the chains.

`Mass` numeric vector with the neutral masses.

searchIsotopes

*Target isotopes search***Description**

This function uses annotation results of an unlabelled sample to search for labelled compounds in a labelled sample.

Usage

```
searchIsotopes(
  msobject,
  label,
  adductsTable = LipidMS::adductsTable,
  rttol = 10,
  ppm = 10,
  coelCutoff = 0.8
)
```

Arguments

<code>msobject</code>	a msobject generated by any of the identification functions.
<code>label</code>	isotope employed for the experiment. It can be "13C" or "D".
<code>adductsTable</code>	adducts table employed for lipids annotation.
<code>rttol</code>	rt window in seconds.
<code>ppm</code>	mass error tolerance.
<code>coelCutoff</code>	coelution score threshold between isotopes. By default, 0.8.

Value

List with the isotopes for each compound in the results data frame.

Author(s)

M Isabel Alcoriza-Balaguer <maribel_alcoriza@iislafe.es>

smdb	<i>SMs database</i>
------	---------------------

Description

In silico generated database for common SMs.

Usage

```
data("smdb")
```

Format

Data frame with 52 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

sphdb	<i>Sphingoid bases database</i>
-------	---------------------------------

Description

In silico generated database for common sphingoid bases.

Usage

```
data("sphdb")
```

Format

Data frame with 4 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

sphPdb

Sphingoid bases phosphate database

Description

In silico generated database for common sphingoid bases phosphate.

Usage

```
data("sphPdb")
```

Format

Data frame with 4 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

tgdb

TGs database

Description

In silico generated database for common TGs.

Usage

```
data("tgdb")
```

Format

Data frame with 376 observations and the following 3 variables.

formula character vector containing molecular formulas.

total character vector indicating the total number of carbons and double bounds of the chains.

Mass numeric vector with the neutral masses.

Index

* datasets

- adductsTable, [3](#)
- baconjdb, [4](#)
- badb, [5](#)
- carnitinesdb, [5](#)
- CEdb, [6](#)
- cerdb, [6](#)
- cldb, [12](#)
- confLevels, [16](#)
- dgdb, [21](#)
- fadb, [21](#)
- fahfadb, [22](#)
- hfadb, [25](#)
- lysopadb, [97](#)
- lysopcdb, [97](#)
- lysopedb, [98](#)
- lysopgdb, [98](#)
- lysopidb, [99](#)
- lysopsdb, [99](#)
- mgdb, [100](#)
- nlsphdb, [100](#)
- padb, [102](#)
- pcdb, [103](#)
- pedb, [103](#)
- pgdb, [104](#)
- pidb, [104](#)
- psdb, [106](#)
- smdb, [107](#)
- sphdb, [107](#)
- sphPdb, [108](#)
- tgdb, [108](#)
- adductsTable, [3](#), [23](#)
- assignDB, [4](#), [17](#), [26](#), [28](#), [30](#), [32](#), [35](#), [37](#), [40](#), [42](#),
[44](#), [46](#), [48](#), [50](#), [52](#), [54](#), [56](#), [58](#), [60](#), [62](#),
[64](#), [66](#), [68](#), [71](#), [73](#), [76](#), [78](#), [80](#), [82](#), [84](#),
[86](#), [88](#), [90](#), [92](#), [95](#)
- baconjdb, [4](#)
- badb, [5](#)
- carnitinesdb, [5](#)
- CEdb, [6](#)
- cerdb, [6](#)
- chainFragments, [7](#), [7](#), [15](#), [26](#), [28](#), [30](#), [32](#), [35](#), [37](#),
[39](#), [42](#), [46](#), [48](#), [50](#), [52](#), [54](#), [56](#), [58](#), [63](#),
[66](#), [68](#), [71](#), [73](#), [75](#), [78](#), [82](#), [84](#), [94](#), [95](#)
- checkClass, [9](#), [28](#), [30](#), [32](#), [34](#), [35](#), [37](#), [39](#), [42](#),
[44](#), [46](#), [48](#), [50](#), [52](#), [54](#), [56](#), [58](#), [60](#), [63](#),
[66](#), [68](#), [70](#), [71](#), [73](#), [75](#), [78](#), [81](#), [82](#), [84](#),
[86](#), [88](#), [90](#), [92](#), [94](#), [101](#)
- checkIntensityRules, [10](#), [32](#), [35](#), [37](#), [39](#), [40](#),
[42](#), [63](#), [64](#), [66](#), [68](#), [71](#), [73](#), [75](#), [76](#), [78](#),
[82](#), [84](#), [95](#), [101](#)
- cldb, [12](#)
- coelutingFragments, [7](#), [9](#), [13](#)
- coelutionScore, [14](#)
- combineChains, [7](#), [11](#), [14](#), [101](#)
- confLevels, [16](#)
- createLipidDB, [7](#), [9](#), [16](#), [17](#), [23](#), [26](#), [28](#), [30](#),
[32](#), [35](#), [37](#), [40](#), [42](#), [44](#), [46](#), [48](#), [50](#), [52](#),
[54](#), [56](#), [58](#), [60](#), [62](#), [64](#), [66](#), [68](#), [71](#), [73](#),
[76](#), [78](#), [80](#), [82](#), [84](#), [86](#), [88](#), [90](#), [92](#), [95](#)
- crossTables, [17](#)
- dataProcessing, [13](#), [18](#), [23](#), [26](#), [28](#), [30](#), [32](#),
[34](#), [37](#), [39](#), [41](#), [44](#), [46](#), [48](#), [50](#), [52](#), [54](#),
[56](#), [58](#), [60](#), [61](#), [63](#), [65](#), [68](#), [70](#), [73](#), [75](#),
[77](#), [80](#), [81](#), [84](#), [86](#), [88](#), [90](#), [92](#), [94](#)
- ddaFragments, [20](#)
- dgdb, [21](#)
- fadb, [21](#)
- fahfadb, [22](#)
- findCandidates, [9](#), [13](#), [15](#), [20](#), [22](#), [23](#), [101](#)
- getInclusionList, [24](#)
- hfadb, [25](#)
- idBANeg, [25](#)
- idCarpos, [27](#)

idCEpos, [29](#)
idCerneg, [31](#)
idCerpos, [34](#)
idCLneg, [36](#)
idDGpos, [38](#)
idFAHFAneg, [41](#)
idFAneg, [43](#)
idLPCneg, [45](#)
idLPCpos, [47](#)
idLPEneg, [49](#)
idLPEpos, [51](#)
idLPGneg, [53](#)
idLPIneg, [55](#)
idLPSneg, [57](#)
idMGpos, [59](#)
idNEG, [61](#)
idPCneg, [62](#)
idPCpos, [65](#)
idPENeg, [67](#)
idPEpos, [70](#)
idPGneg, [72](#)
idPGpos, [74](#)
idPINeg, [77](#)
idPOS, [79](#)
idPSneg, [81](#)
idSMpos, [83](#)
idSphneg, [85](#)
idSphPneg, [87](#)
idSphpos, [89](#)
idSphPpos, [91](#)
idTGpos, [93](#)

LipidMSapp, [96](#)
lysopadb, [97](#)
lysopcdb, [97](#)
lysopedb, [98](#)
lysopgdb, [98](#)
lysopidb, [99](#)
lysopsdb, [99](#)

mgdb, [100](#)

nlsphdb, [100](#)

organizeResults, [101](#)

padb, [102](#)
pcdb, [103](#)
pedb, [103](#)

pgdb, [104](#)
pidb, [104](#)
plotLipids, [105](#)
psdb, [106](#)

searchIsotopes, [106](#)
smdb, [107](#)
sphdb, [107](#)
sphPdb, [108](#)

tgdb, [108](#)